

# MedtronicCSVtoFHIR

MedtronicCSVtoFHIR is a Python SCRIPT that facilitates the conversion of Medtronic CSV data to FHIR (Fast Healthcare Interoperability Resources) bundles for interoperability with healthcare systems.

## Description

This library processes Medtronic CSV data containing patient information and generates FHIR-compliant data bundles. It includes functions for handling glucose observations, carbohydrate observations, and insulin medication administration, allowing seamless integration of data into healthcare systems.

## Features

- Conversion of Medtronic CSV data to FHIR bundles ( in [main.py](#)).
- Generation of glucose observations ( in [conversion.py](#)).
- Generation of carbohydrate observations ( in [conversion.py](#)).
- Creation of insulin medication administration bundles ( in [conversion.py](#)).
- Integration with FHIR-compatible systems [VERSION R4] (customizable via environment variables (at [.env](#)))

## DEMO

We've included anonymized patient data produced between October and November 2023 in a separate folder for testing purposes.

The data is in the [Demo folder](#).

Additionally, associated FHIR bundles have been generated based on the provided data using the following format:

- Glucose bundles (e.g., [bundleyear2023month10week1part1.json](#), )
- Carbohydrate bundles (e.g., [carbsbundleyear2023month10week1part\\_1](#) )
- Medication Administration bundles (Insulin) (e.g., [insulinbundleyear2023month11week2part\\_2.json](#) )

Note: Some weekly data may be divided into multiple parts, especially larger bundles that could potentially exceed the processing capacity of an FHIR server.

## Setup

1. Clone the repository:

```
```bash git clone https://github.com/pr008/Medtronic-csv-FHIR
cd Medtronic-csv-FHIR ```
```

1. Run the installation script for macOS (and other Unix-like systems like iOS):

```
bash ./install.sh
```

1. Run the installation script for Windows:

```
bash install.bat
```

This script will perform the following tasks:

- Checks if Python is installed. If not, it provides instructions to download and install Python.
- Creates a virtual environment named venv in the project directory.
- Activates the virtual environment.
- Installs the required dependencies specified in requirements.txt.
- Once the setup is complete, you can start using the project by following the usage instructions or running your desired commands.

## Installation Step by Step

### Python Installation (if not installed)

If you haven't installed Python, you can download it from the official [a https://www.python.org/downloads/](https://www.python.org/downloads/) and follow the installation instructions based on your operating system.

### Setting up a Virtual Environment (optional but recommended)

While not mandatory, it's recommended to use a virtual environment to isolate project dependencies. To set up a virtual environment:

1. Open a terminal or command prompt.
2. Navigate to the project directory.

```
bash cd /path/to/Medtronic-csv-FHIR
```

1. Run the following command to create a virtual environment (replace venv with your preferred environment name):

```
bash python3 -m venv venv
```

1. Activate the virtual environment:
  - On Windows:

```
bash venv\Scripts\activate
```

- On macOS and Linux:

```
bash source venv/bin/activate
```

## Installing Dependencies

Once Python is installed and a virtual environment (if desired) is set up, you can install the required packages:

1. Ensure you're in the project directory and the virtual environment is activated.
2. Run the following command to install the required dependencies using pip:

```
bash pip install -r requirements.txt
```

This will install all the necessary dependencies specified in the requirements.txt file for the Medtronic CSV to FHIR conversion library.

## Environment Variables and Customization

This will install all the necessary dependencies specified in the requirements.txt file for the Medtronic CSV to FHIR conversion library.

## Environment Variables and Customization

Please be aware that some of these constants are defined by Medtronic and they **MAY CHANGE**

Refer to the [Medtronic CSV Documentation](#) for more details.

Please note that the documentation for Medtronic CSV files is incomplete. As a result, some variables were guessed based on the available information.

For example, for glucose data: - BGSENTFOR\_CALIB: They use these values and report the manual glucose values

- USERACCEPTEDREMOTE\_BG: They use these values and report the manual glucose values

## Medtronic specific code in CSV

**Note:** The official documentation lacks explicit details for variables that distinguish between BOLUS, BASAL, and CORRECTION types in Medtronic data. We've defined several constants to fill this gap. These constants are crucial for interpreting the data accurately, particularly in the auto\_bolus column and similar fields. Here's how we've defined and used these variables:

Variable Name	Description
BG_SENT_FOR_CALIB	Represents manual glucose values, used where <code>df['BG Source'] == 'BG_SENT_FOR_CALIB'</code> .    USER_ACCEPTED_REMOTE_BG
	Indicates remote BG acceptance, used where <code>df['BG Source'] == 'USER_ACCEPTED_REMOTE_BG'</code> .    MEDTRONIC_BG_READIN_RECEIVED
	Corresponds to BG readings received, used where <code>df['BG Source'] == 'BG_READIN_RECEIVED'</code> .    MEDTRONIC_CLOSED_LOOP_AUTO_BASAL
	Part of BASAL_CONSTANT, indicates automated basal insulin delivery.    MEDTRONIC_CLOSED_LOOP_AUTO_INSULIN
	Also part of BASAL_CONSTANT, signifies automated insulin delivery not specified as basal or bolus.    MEDTRONIC_CLOSED_LOOP_BG_CORRECTION_AND_FOOD_BOLUS
	Part of BOLUS_CONSTANT, used for food bolus and BG correction combined.    MEDTRONIC_CLOSED_LOOP_AUTO_BOLUS
	Part of CORRECTION_CONSTANT, used for automated bolus corrections.    MEDTRONIC_CLOSED_LOOP_BG_CORRECTION
	Also part of CORRECTION_CONSTANT, indicates a BG correction without a food bolus.    MEDTRONIC_BOLUS_WIZARD
	Represents bolus calculations or suggestions provided by the Bolus Wizard feature.

These variables are set as environment variables in our code, allowing for dynamic adjustment and easier management of the constants used in data processing.

## Glucose Scale and Codes for FHIR Resources

Variable Name	Description	Usage
GLUCOSE_SCALE_MMOL	Measurement unit: mmol/L	Measurement unit for glucose in mmol/L.    GLUCOSE_SCALE_MMOL_SYSTEM
	System URL for mmol/L scale	URL of the system for mmol/L glucose scale.    GLUCOSE_SCALE_MMOL_CODE
	Code for mmol/L glucose scale	Code for mmol/L glucose scale.    GLUCOSE_SCALE_MMOL_DISPLAY
	Display name for mmol/L glucose scale	Display name for mmol/L glucose scale.    GLUCOSE_SCALE_MG
	Measurement unit: mg/dL	Measurement unit for glucose in mg/dL.    GLUCOSE_SCALE_MG_SYSTEM
	System URL for mg/dL scale	URL of the system for mg/dL glucose scale.    GLUCOSE_SCALE_MG_CODE
	Code for mg/dL glucose scale	Code for mg/dL glucose scale.    GLUCOSE_SCALE_MG_DISPLAY
	Display name for mg/dL glucose scale	Display name for mg/dL glucose scale.    CGM_GLUCOSE_SYSTEM
	System URL for CGM glucose	URL of the system for CGM glucose.    CGM_GLUCOSE_CODE
	Code for CGM glucose measurement	Code for CGM glucose measurement.    CGM_GLUCOSE_DISPLAY
	Display name for CGM glucose measurement	Display name for CGM glucose measurement.

## Manual Blood Glucose and Interpretation for FHIR Resources

Variable Name	Description	Usage
BG_GLUCOSE_SYSTEM	System URL for manual blood glucose	URL of the

system for manual blood glucose. | | BG\_GLUCOSE\_CODE | Code for manual blood glucose | Code for manual blood glucose. | | BG\_GLUCOSE\_DISPLAY | Display name for manual blood glucose | Display name for manual blood glucose. | | GLUCOSE\_INTERPRETATION\_SYSTEM | System URL for glucose interpretation | URL of the system for glucose interpretation. | | GLUCOSE\_INTERPRETATION\_CODE | Code for glucose interpretation | Code for glucose interpretation. | | GLUCOSE\_INTERPRETATION\_DISPLAY | Display name for glucose interpretation | Display name for glucose interpretation. | | GLUCOSE\_INTERPRETATION\_LU\_CODE | Code for significantly low glucose | Code for significantly low glucose interpretation. | | GLUCOSE\_INTERPRETATION\_LU\_DISPLAY | Display name for significantly low glucose | Display name for significantly low glucose. | | GLUCOSE\_INTERPRETATION\_LU\_MMOL | Significantly low glucose value (mmol/L) | Significantly low glucose value in mmol/L. | | GLUCOSE\_INTERPRETATION\_LU\_MG | Significantly low glucose value (mg/dL) | Significantly low glucose value in mg/dL. | | GLUCOSE\_INTERPRETATION\_L\_CODE | Code for low glucose | Code for low glucose interpretation. | | GLUCOSE\_INTERPRETATION\_L\_DISPLAY | Display name for low glucose | Display name for low glucose interpretation. | | GLUCOSE\_INTERPRETATION\_L\_MMOL | Low glucose value (mmol/L) | Low glucose value in mmol/L. | | GLUCOSE\_INTERPRETATION\_L\_MG | Low glucose value (mg/dL) | Low glucose value in mg/dL. | | GLUCOSE\_INTERPRETATION\_N\_CODE | Code for normal glucose | Code for normal glucose interpretation. | | GLUCOSE\_INTERPRETATION\_N\_DISPLAY | Display name for normal glucose | Display name for normal glucose interpretation. | | GLUCOSE\_INTERPRETATION\_H\_CODE | Code for high glucose | Code for high glucose interpretation. | | GLUCOSE\_INTERPRETATION\_H\_DISPLAY | Display name for high glucose | Display name for high glucose interpretation. | | GLUCOSE\_INTERPRETATION\_H\_MMOL | High glucose value (mmol/L) | High glucose value in mmol/L. | | GLUCOSE\_INTERPRETATION\_H\_MG | High glucose value (mg/dL) | High glucose value in mg/dL. | | GLUCOSE\_INTERPRETATION\_HU\_CODE | Code for significantly high glucose | Code for significantly high glucose interpretation. | | GLUCOSE\_INTERPRETATION\_HU\_DISPLAY | Display name for significantly high glucose | Display name for significantly high glucose. | | GLUCOSE\_INTERPRETATION\_HU\_MMOL | Significantly high glucose value (mmol/L) | Significantly high glucose value in mmol/L. | | GLUCOSE\_INTERPRETATION\_HU\_MG | Significantly high glucose value (mg/dL) | Significantly high glucose value in mg/dL. |

## Insulin Codes (Bolus and Basal) for FHIR Resources

You can customize these as well.

Variable Name	Description	Usage
MEDICATION_ADMINISTRATION_SYSTEM	System URL for bolus insulin system	System URL for background insulin (bolus).
MEDICATION_ADMINISTRATION_CODE		
MEDICATION_ADMINISTRATION_DISPLAY		

MEDICATION\_ADMINISTRATION\_BOLUS\_CODE | Code for background insulin (bolus) | Code for background insulin (bolus). | |  
 MEDICATION\_ADMINISTRATION\_BOLUS\_DISPLAY | Display name for background insulin | Display name for background insulin. | |  
 MEDICATION\_ADMINISTRATION\_BOLUS\_UNIT\_SYSTEM | URL for bolus insulin unit system | System URL for insulin unit (bolus). | |  
 MEDICATION\_ADMINISTRATION\_BOLUS\_UNIT\_CODE | Code for insulin unit (bolus) | Code for insulin unit (bolus). | |  
 MEDICATION\_ADMINISTRATION\_BOLUS\_UNIT\_DISPLAY | Display name for insulin unit (bolus) | Display name for insulin unit (bolus). | |  
 MEDICATION\_ADMINISTRATION\_BASAL\_SYSTEM | URL for basal insulin system | System URL for bolus insulin. | | MEDICATION\_ADMINISTRATION\_BASAL\_CODE | Code for bolus insulin | Code for bolus insulin. | |  
 MEDICATION\_ADMINISTRATION\_BASAL\_DISPLAY | Display name for bolus insulin | Display name for bolus insulin. | |  
 MEDICATION\_ADMINISTRATION\_BASAL\_UNIT\_SYSTEM | URL for basal insulin unit system | System URL for insulin unit (basal). | |  
 MEDICATION\_ADMINISTRATION\_BASAL\_UNIT\_CODE | Code for insulin unit (basal) | Code for insulin unit (basal). | |  
 MEDICATION\_ADMINISTRATION\_BASAL\_UNIT\_DISPLAY | Display name for insulin unit (basal) | Display name for insulin unit (basal). |

## Carbohydrates (Estimated)

Variable Name	Description	Usage
CARBOHYDRATES_EST_SYSTEM	URL for carbohydrate intake system	System URL for estimated carbohydrate intake.
CARBOHYDRATES_EST_CODE	Code for carbohydrate intake	Code for estimated carbohydrate intake.
CARBOHYDRATES_EST_DISPLAY	Display name for carbohydrate intake	Display name for estimated carbohydrate intake.
CARBOHYDRATES_EST_UNIT	Unit of measurement for carbohydrates	Unit of measurement for estimated carbohydrates.
CARBOHYDRATES_EST_UNIT_SYSTEM	URL for carbohydrate unit system	System URL for carbohydrate unit.
CARBOHYDRATES_EST_UNIT_CODE	Code for carbohydrate unit	Code for estimated carbohydrate unit.